

PRODUCT RELEASE CERTIFICATE

DECON-CLEAN

Date of Manufacture: 17-Oct-2025

Product DECON-CLEAN

Fill Size Gallon Simple Mix

Cat No: DC-04-1Z

Lot #: 25-DCSM-522179L

Expiration: 01-Oct-2027

TEST	SPECIFICATION	RESULT
pH of Concentrate	9.0 - 11.0	10.0
Solubility in Water	Completely	PASS
Specific Gravity (@20°C)	1.010 - 1.040	1.014
Appearance	Clear/Free of Suspended Matter	PASS
% Non-Volatile Matter	≤ 20.0%	10.7%
VAI WFI Water		
Nitrates	< 0.2 ppm	< 0.2 ppm
Conductivity	≤ 5 uS/cm	1.0 uS/cm
TOC	< 8.0 mg/L	0.054 mg/L
Ammonia	≤ 0.2 ppm	< 0.2 ppm
Appearance	Clear Liquid	PASS
Residue on Evaporation	Max 0.001%	0.000%
pH	5.0 - 7.0	5.9
Calcium	No Turbidity	PASS
LAL Result	< 0.25 EU/mL	< 0.015 EU/mL
Sterility		
USP Sterility	PASS	PASS
Gamma Irradiation		
Gamma Irradiation Dose Min	25.0 - 50.0 kGy	26.5 kGy
Gamma Irradiation Dose Max	25.0 - 50.0 kGy	46.7 kGy



Authorized Signature - Quality Assurance

20 Nov 2025

Date

CONTAMINATION CONTROL INNOVATIONS SINCE 1981

Corporate Headquarters 15 Lee Boulevard, Malvern, PA 19355-1234 USA, T: 610-644-8335 F: 610-644-8336
www.sterile.com