



Model:
LED200 - 631-1085

Typology:
ILLUMINATOR

Description:
P-LED cold light generator

Light source	High power 14W P-LED
Chromaticity coordinates (x,y)	0.320, 0.320 (pure white, 6300K)
Optical power (with double-arm)	Illuminance (at 10 cm distance): 70000 lux
Power supply	Brightness control potentiometer. External power supply provided. Input: 110/240 Vac 50/60 Hz; output: 6 Vdc 5A
Dimensions	Length: 200 mm Width: 120 mm Weight: 4,1 Kg
Standard configuration:	Double-arm optical fiber guide. Each arm fitted with focusable lens adapter. Arm length 500 mm.