



Manufacturer's Declaration of Conformity
(in accordance with ISO/IEC 17050-1)

Issuer's name: Tuttnauer Europe B.V.
Issuer's address: Hoeksteen 11
P.O. Box 7191
4800 GD Breda
The Netherlands

Object of the declaration: Laboratory Autoclaves
Model Codes: See Addendum: product list

In addition we ensure and declare that the distributed products, as mentioned, meet the relevant provisions.

Document	Title
- Directive 2014/35/EU	Low Voltage Directive
- Directive 2014/30/EU	EMC Directive
- EN 61326-1:2013	Electrical equipment for measurement, control and laboratory use – EMC requirements
- EN 61010-1:2010	Safety requirements for electrical equipment for measurement, control and laboratory use
- EN 61010-2-040:2015	Particular requirements for sterilizers and washer-disinfectors used to treat medical materials

All harmonized standards which are applicable to the object, as published in the Official Journal of the European Communities.

Additional information:

Notified Body: Dekra Certification BV
Meander 1051
6825 MJ Arnhem
The Netherlands

First issue date certificate: February 14, 2010
Valid until: February 1, 2021

This documentation is supported by the Quality Management System certification based on the harmonized standard ISO 9001:2015, Quality Management System Certificate with reference numbers: 2129743, issued for the 1st time on February 14, 2010, and delivered by Dekra Certification BV.

Signed for and on behalf of: Breda, The Netherlands, February 06, 2018
Name: A. Geilings
Function (Company): Managing Director

Signature:



ADDENDUM: Product List

Laboratory Series

Model	Inner Dimensions Diameter x Depth (mm)	Volume Liters
EL-D Series - Fully Automatic		
2840EL-D	280 x 400	28
3840EL-D	380 x 400	52
3850EL-D	380 x 500	65
3870EL-D	380 x 690	85
5050EL-D	500 x 500	110
5075EL-D	500 x 750	160
ELC-D Series - Fully Automatic with Rapid Cooling		
3840ELC-D	380 x 400	52
3850ELC-D	380 x 500	65
3870ELC-D	380 x 690	85
5050ELC-D	500 x 500	110
5075ELC-D	500 x 750	160
ELV-D Series - Fully Automatic, Vertical		
2840ELV-D	280 x 460	31
3840ELV-D	380 x 400	52
3850ELV-D	380 x 500	65
3870ELV-D	380 x 690	85
5050ELV-D	500 x 500	110
5075ELV-D	500 x 750	160
ELVC-D Series - Fully Automatic, Vertical, with Rapid Cooling		
2840ELVC-D	280 x 460	31
3840ELVC-D	380 x 400	52
3850ELVC-D	380 x 500	65
3870ELVC-D	380 x 690	85
5050ELVC-D	500 x 500	110
5075ELVC-D	500 x 750	160