

PRODUCT SPECIFICATION SHEET

Spec. NEU - 01/01/2014

Product Name:	Decon NEUTRACON In certain countries, Decon Neutracon is sold under the trade name NEUTRAD.
Description:	Surface active cleaning agent, for laboratory, medical and industrial applications. NOT suitable for <i>domestic</i> applications. Supplied as a liquid concentrate, for dilution with water. For 'manual' cleaning, or for use in an ultrasonic cleaning bath.
Properties:	Biodegradable*, phosphate-free, non-toxic, non-flammable, non-corrosive and totally rinsable.
Appearance:	Pale amber viscous liquid, having faint odour.
pH:	Decon NEUTRACON is a 'near neutral' concentrate. pH 7.5 - 1% solution in water (typical).
Packaging:	Supplied in <i>1 litre, 5 litre and 20 litre</i> high density polyethylene containers, with tamper-evident screw-on caps.
Applications:	For cleaning and/or decontaminating non-ferrous metals, such as aluminium, zinc, copper, brass, etc., 'soft' glass, coated lenses, polymer surfaces and other sensitive surfaces which must not be cleaned by acidic or alkaline cleaning agents.
Formulation:	A mixture of anionic and non-ionic surface active agents, in an aqueous polyhydric alcohol base.

*** Ready Biodegradability by OECD 301E (ISO method 7827 - 1986 (E))**

© Copyright: Decon Laboratories Limited U.K. : 1997



Decon Laboratories Limited
Conway Street, Hove, East Sussex, BN3 3LY, England

