



Sol - Care® Safety Blood Collection Needle

- Transparent design that allows visualization of the blood flashback and flow
- Single-handed activation option that provides more flexibility to care for the patient and venipuncture site
- Sliding safety shield pushes over the needle and locks with an audible click, protecting against needlestick injury

***samples available upon request**

VWR Cat No.

Product Description

76582-466 Sol-Care Safety Blood Collection Needle 21G x 3/4" x 7"

76582-468 Sol-Care Safety Blood Collection Needle 21G x 3/4" x 12"

76582-470 Sol-Care Safety Blood Collection Needle 23G x 3/4" x 7"

76582-472 Sol-Care Safety Blood Collection Needle 23G x 3/4" x 12"

76582-474 Sol-Care Safety Blood Collection Needle 25G x 3/4" x 7"

76582-476 Sol-Care Safety Blood Collection Needle 25G x 3/4" x 12"

VWR Cat No.

Product Description

76582-456 Sol-Care Safety Blood Collection Needle 21G x 3/4" x 7" w/ pre-attached Tube Holder

76582-454 Sol-Care Safety Blood Collection Needle 21G x 3/4" x 12" w/ pre-attached Tube Holder

76582-460 Sol-Care Safety Blood Collection Needle 23G x 3/4" x 7" w/ pre-attached Tube Holder

76582-458 Sol-Care Safety Blood Collection Needle 23G x 3/4" x 12" w/ pre-attached Tube Holder

76582-464 Sol-Care Safety Blood Collection Needle 25G x 3/4" x 7" w/ pre-attached Tube Holder

76582-462 Sol-Care Safety Blood Collection Needle 25G x 3/4" x 12" w/ pre-attached Tube Holder