

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Sample Library Tubes and Closures, 1.2 ml
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2025-12-22

VWR Catalogue Number	Description	Unit
82006-698	Individual Tubes in 8 x 12 Racks, 1.2 ml	Pack of 960
		Case of 9600
93000-040	Individual Tubes in 8 x 12 Racks, Sterile, 1.2 ml	Pack of 960
		Case of 9600
82006-700	8-Strip Tubes in 8 x 12 Racks, 1.2 ml	Pack of 960
		Case of 9600
82006-702	12-Strip Tubes in 8 x 12 Racks, 1.2 ml	Pack of 960
		Case of 9600
93000-044	12-Strip Tubes in 8 x 12 Racks, Sterile, 1.2 ml	Pack of 960
		Case of 9600
83009-678	Individual Tubes, 1.2 ml	Pack of 1,000
		Case of 10000
83009-680	8-Strip Tubes, 1.2 ml	Pack of 125
		Case of 1250
83009-683	12-Strip Tubes, 1.2 ml	Pack of 125
		Case of 1250

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with the **ISO 9001:2015** standard.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Polypropylene (PP)
Sterilization	The products are non-sterile.
ATP Assay	The products are free from ATP.
DNase & RNase Free	The products are free of any detectable DNase/RNase contamination. Test sensitivity: Deoxyribonuclease (DNase): < 10 ⁻⁷ Kunitz Units/μl Ribonuclease (RNase): < 10 ⁻⁹ Kunitz Units/μl. Human and mouse DNA: <30pg
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin, and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Based on the resin manufacturer the products are free from toxins.
Non-Pyrogenic Statement	Not tested.
Latex Statement	Based on the resin manufacturer the products are free from latex.
BPA Statement	Based on the resin manufacturer the products are free from Bisphenols.
DEHP Statement	Based on the resin manufacturer the products are free from DEHP.
RoHS	Not Applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routine analysis is performed.
Storage Conditions	Storage in typical warehouse or laboratory temperatures is suitable.
Shelf-Life	5 years.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	December 2025