

03/04/2025

Freudenberg Medical Europe
2-8 Leibigstrasse
67661 KAISERSLAUTERN
GERMANY

PO# 4500039832

QTY. 4 EA

PRODUCT CERTIFICATION


<u>FREUDENBERG P/N</u>	<u>DESCRIPTION</u>	<u>FREUDENBERG L/N</u>
70-400-07		0000329994

**¾" SILICONE GASKET
DOW CORNING, SILASTIC Q7-4780 MATERIAL**

<u>MATERIAL</u>	<u>FM P/N</u>	<u>FM L/N</u>
Silicone	30-011-01	0000317092
Elastomer	30-012-01	0000317093

Freudenberg Medical certifies that our silicone gaskets are manufactured following Good Manufacturing Practices (GMP) in accordance with the new Quality System Regulation (QSR) in a FDA registered medical device facility. Freudenberg Medical has been certified to meet ISO 13485 quality system. Material manufacturer states that these materials meet or exceed USP Class VI biocompatibility requirements and FDA Regulations CFR 177.2600 for contact with food and drugs. These platinum-cured silicone products have full traceability to raw materials and contain no peroxide catalysts.

Sincerely,


Franco Castrejon / Shipping Clerk

(7040007 REV. B)

Freudenberg Medical
Material Number Batch Number

45639002-0000329994