

Material Number: 87039.260
Material description: Ion chromatography anion multi-element standard solution I
Grade: ARISTAR®
Date specification set: 25th Sep 2019

(All batches of this product after this date are tested to this specification. In the unlikely event that the specification is changed a new certificate will be issued)

Characteristics	Specifications	Traceability
F ⁻	99.8 – 100.2 mg/l	NIST SRM No 3183
Cl ⁻	249.5 – 250.5 mg/l	NIST SRM No 3182
NO ₃ ⁻	499 – 501 mg/l	NIST SRM No 3185
SO ₄ ²⁻	499 – 501 mg/l	NIST SRM No 3181
PO ₄ ³⁻	998 – 1002 mg/l	NIST SRM No 3186

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International