

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Flasks, Flat Bottom, with Standard Ground Joint
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2025-03-14

Catalogue Number	Capacity [mL]	Height [mm]	Ext.Ø [mm]	Socket size	Pk
201-1344	50	85	51	24/29 NS	2
201-1380	50	85	51	29/32 NS	2
201-1348	100	103	64	24/29 NS	2
201-1381	100	103	64	29/32 NS	2
201-1360	250	130	85	24/29 NS	2
201-1382	250	130	85	29/32 NS	2
201-1361	500	160	105	24/29 NS	2
201-1383	500	160	105	29/32 NS	2
201-1384	1000	187	131	29/32 NS	1

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 17025:2017** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirements of the following standards and regulations

Technical standards for products: DIN 12353

Calibration procedure: NA

Material	The products are made from DIN ISO 3585:1998, TYPE-1, CLASS- A, Boro 3.3 Glass that meet ASTM E-438, TYPE-1, CLASS- A Standard.
Sterilization	Not applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable
Storage Conditions	Store at room temperature in normal warehouse conditions.
Shelf Life	No predefined shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	03-2025