

Material	26843.268
Material description	Potassium iodide
Grade	GPR RECTAPUR
CAS Number	7681-11-0
Molecular formula	KI
Molecular mass	166.00

Characteristics	Specifications
Assay	$\geq 99 \%$
Heavy metals (as Pb)	$\leq 10 \text{ ppm}$
Cl (Chloride)	$\leq 0.1 \%$
Fe (Iron)	$\leq 10 \text{ ppm}$

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.