




VWR® Frosted Microscope Slides without Tissue Paper

VWR Frosted Microscope Slides without Tissue Paper are specially designed with consistent and clean surfaces for staining experiments. They are ideal for histology, hematology, and microbiology applications. Slides feature frosted areas on both sides of one end for convenient marking. A fine polished edge and 45° clipped corners promote ease and safety in handling. Slides are packaged 72 per pack, 20 packs per case.

Designed for routine experiments where no extra adhesion or no adaptability to auto-writers is required for blood or body liquid tests, frozen tissue section examinations, biological specimen preparations, etc., in histology-pathology, hematology, cytology, microbiology, and other similar applications.

- Non-Sterile
- Store in a cool, dry place
- Conforms with the essential requirements of European in vitro Diagnostic 98/ 79 EC.
- Term of Validity: 18 months

Specifications:

Color	Transparent
Material	Green glass
Component	SiO ₂ (~70-73%):Na ₂ O(~13-15%):CaO(~7-12%):MgO(~1-4.5%):Al ₂ O ₃ (~1-2%):SO ₃ (~0-0.3):Fe ₂ O ₃ (~0.05-0.2%)
Dimensions	75.0(+0.5) x 25.0(+0.5)
Edge	Ground edges
Corner	45°
Marking Area	Double sides frosted at one end
Tissue Paper	Without tissue paper
Surface Modification	No modification
Thickness (mm)	1.1±0.1
Gross Weight (kg)	8.5
Net Weight (kg)	7.5
Refractivity	1.50-1.54

VWR® Frosted Microscope Slide

Size	Thickness	Pack of	Case of	Cat. No.
75 X 25 mm	1 mm	72	20 packs	89085-399