

Product Data Sheet for KIMTECH SCIENCE* Precision Wipes - White /Small - 7552

Product Information

Product Long Code	07552000
Customer Facing Code	7552
Description	KIMTECH SCIENCE* Precision Wipes - White /Small
PDS Reference	21/10/2020



Packaging Configuration

Case Contents	7552: 30 Dispenser Boxes x 280 Sheets = 8400 Sheets
Outer Pack Material	Corrugated

Icon	Unit	Dimension					Weight(kg)	
	Dispenser Box	12.00	x	7.30	x	11.80	Length x Width x Height(cm)	
	Sheet	21.30	x	11.40			Length x Width(cm)	

Product Description

Effective cleaning of equipment, materials and working environments is important to ensure the accuracy, stability and longevity of your lab. Kimtech™ SCIENCE provides wipes for delicate tasks, designed to maintain cleanliness in your laboratory or working environment. The light duty, 1-ply wipes are suitable for a variety of precision wiping tasks such as the cleaning of delicate glassware, equipment and surfaces. Lenses, slides and optical instruments can be wiped with ease, and the wipes are suitable for ISO Class 4 or higher laboratory environments. Made from 100% pure fibre, giving high absorbency and a soft feel, the white-coloured wipes are chemically inert allowing use with solvents and chemicals. The wipes are also dispensed from a flexible pop-up box, to ensure clean and fast dispensing, even while on the move. Suitable for use within laboratories in manufacturing and healthcare, as well as industry sectors including optical, laser and electronics, the wipes feature low-lint tissue softness, making them ideal for delicate cleaning applications. They are also suitable for cleaning liquids and dust, and are equipped with anti-static dispensing to reduce lint and antistatic discharge. 30 dispenser boxes x 280 white, small 1 ply sheets = 8400 sheets

Product Features

- 30 dispenser boxes x 280 white, small 1 ply sheets = 8400 sheets
- For delicate equipment, surfaces, optics and materials, our low-lint and high absorbency white wipes have you covered.
- Light duty, white-coloured wipes made from 100% pure fibre
- Provided in a hygienic, user-friendly pop-up box
- Low-lint absorbency for delicate applications and surfaces
- Suitable for use in ISO Class 4 or higher environments
- Suitable for applications in healthcare, manufacturing, optics and electronics. Chemically inert and anti-static

Product Specifications

Fold Type	I	Impregnated	No
Bar Code (Case)	00036000755206	Bar Code (Pack)	0360005515
Recycled Fibre(%)	3	Virgin Fibre(%)	97
Bleaching Method	Elemental Chlorine Free (ECF)	Colour	White

Product Data Sheet for KIMTECH SCIENCE* Precision Wipes - White /Small - 7552

International Standards and Symbols

Solution

Safer Workplace

We know staff safety is your top priority. You can rely on our solutions to protect them, their environment and your important work.



Productive Workplace

The success of your business hinges on your work efficiency. That's why our solutions are designed to help improve your productivity.



Packaging

Unit Of Sale	Case	Net Weight(kg)	3.9
Max Weight(kg)	5.85	Case Dimension (LXWXH)(cm)	36.9 x 24.8 x 37.9
Inner Pack Material	Cartonboard		

Pallet Type	LxWxH (m)	Cases per Layer	Layers per Pallet	Cases per Pallet
S2	1.20 x 1.00 x 1.39	6	5	30

Batch Code



Information

It is the employer's responsibility to assess the risk of the task to be undertaken and determine the correct choice of wipers for the task. The manufacturer, Kimberly-Clark, does not accept any responsibility for the incorrect choice or misuse of wipers shown in this brochure. All care has been taken to ensure that the information contained herein is as accurate as possible at the time of publication, however errors may occur and legislation concerning personal protective equipment is under constant review and may change in the lifetime of this brochure. Accordingly, the specifications for the products may be subject to change. We would advise you to contact INFOFAX if you have any queries concerning the products shown or the suitability of such products for the particular task. Always dispose of used protective equipment in a safe and appropriate manner in accordance with European, National and Local environmental regulations

Websites www.kcprofessional.co.uk

Customer Service Number For more information please contact Customer Service during regular office hours-0870 551 044

Country of Origin UNITED STATES