



Certificate of Conformity

VWR Drosophila Vial

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Ø ext.×H (mm)	Packed	Pk
734-2255	25×95	Bulk	500
734-2259	25×95	Tray	500
734-2262	28.5×95	Bulk	500
734-2265	28.5×95	Tray	500

Material: Polystyrene (PS): This resin has glass like clarity, allowing very clear, distortion free viewing of fly gender and behaviour.

Storage: Store in a clean, dry environment that is at room temperature

Other Characteristics:

- For storage, handling and transportation of Drosophila flies.
- Manufactured in ISO 9001:2000 certified facility.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **April 2019**