



Item Number HY104
Item Hydrochloric Acid, 37%, NF
CAS Number 7647-01-0
Molecular Formula HCl
Molecular Weight 36.46
Assay (w/w)..... 36.5-38.0%

MAXIMUM LIMITS

Residue on Ignition..... To pass test (limit about 0.008%)
Bromide or Iodide..... To pass test
Free Bromine or Chlorine..... To pass test
Sulfate (SO₄)..... To pass test
Sulfite (SO₃)..... To pass test
Heavy Metals (as Pb)..... 5 ppm
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

<i>Corporate Headquarters:</i>	<i>East Coast Facility:</i>
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	