

Certificate of Conformity

VWR Catalogue Number	710-0442
Description	VWR® M40, Transport Swab
Country of Origin	Manufactured in Italy
Date of Issue (yyyy-mm-dd)	2026-04-17

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 13485:2016**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirement

Product Specifications

Medium Visual Appearance and Mechanical Controls: **passed according to internal SOP.**

pH of the medium at 25 °C acceptability limits 7.2 ± 0.3 : **passed**

Microscopic Bioburden Specifications, maximum 2 micro-organisms per 20 microscopic fields at 10X100 magnification: **passed**.

Cultural Response: swabs were inoculated with the tester organisms (see the table) and placed in the medium. Recovery obtained after 24 hours at 20-25°C: **passed**.

TRANSYSTEM M40, AMIES AND STUART CLEAR GEL AND WITH CHARCOAL medium	TRANSYSTEM M40, AMIES AND STUART LIQUID medium
Neisseria gonorrhoeae ATCC 43069	Neisseria gonorrhoeae ATCC 43069
Haemophilus influenzae ATCC 10211	Haemophilus influenzae ATCC 10211
Streptococcus pyogenes ATCC 19615	Streptococcus pyogenes ATCC 19615
Streptococcus pneumoniae ATCC 6305	Streptococcus pneumoniae ATCC 6305
Bacteroides fragilis ATCC 25285	/
Peptostreptococcus anaerobius ATCC 27337	/
Propionibacterium acnes ATCC 6919	/

Material	Cap: High Density Polyethylene (HDPE). Tube: Polypropylene (PP). Stick Material: Polystyrene (PS). Tip: Rayon.
Sterilization	Products are gamma irradiated. Sterilisation procedure is validated in conformity with BS EN ISO 11137-1:2015 and BS EN ISO 11137-2:2015 to guarantee SAL 10 ⁻⁶ .
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	DEHP has not been intentionally incorporated in the manufacturing of the products.
RoHS	Not applicable.
REACH Statement	Product meets the criteria for authorization under the European Union's REACH regulation (EC) No. 1907/2006.

Storage Conditions 5 – 25 °C.do not freeze.

Shelf Life 15 months.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	05-2026