

Simple and Affordable



Easy Redox Titrations

Everything you need for inflection-point Redox titrations. This versatile titrator has a platinum ring redox-sensing electrode making it the ideal instrument to meet the needs of food, beverage, galvanic and even environmental analyses.



Easy App Guidance

This titrator has been especially designed for simplicity. New operators are already familiar with the Smartphone like Apps oriented user interface. Start your routine analysis with a LongClick™ on the home screen. The intuitive navigation speaks your language and 14 others. Learning and using a new analytical instrument has never been so easy.



Swiss Engineered

EasyPlus™ Titrators are designed and constructed by our Swiss development team according to the highest possible software and hardware standards. Their renowned expertise has made the impossible possible; to develop such a high quality instrument with such a small footprint.



iTitrate for Easy Programming

iTitrate Intelligence makes method programming simple with EasySetup guidance. Set only the minimal parameters and the titrator programs the rest, including correct calculations based on desired measurement units. Easy setup equals swift implementation.



Easy Ox Titrator The Redox Titrator

The Easy Ox Titrator speeds up and simplifies automated reduction and oxidation (redox) titrations. The simple 'App' style user interface reduces training time while increasing work speed and simplicity. iTitrate Intelligence's guided set-up assists in method programming to ensure samples are tested consistently, regardless of operator, and with minimal effort. Eliminate manual calculations and documentation with a printer or EasyDirect™ Titration PC Software.

The Easy Ox Titrator's platinum ring electrode makes it ideal for the following applications and many more:

- Vitamin C determination in food and beverages.
- Peroxide value of cooking oils.
- Reducing sugar content in wines and juices.
- Hydrogen peroxide content of cleaning and plating solutions.

Standard delivery

Product	Quantity	Order Number
EM40-BNC Sensor Redox mV	1	30043106
EasyPlus burette 20 mL	1	30043901
EasyPlus tubing set instrument	1	30065464
EasyStir GT	1	30065467
EasyPlus titration head GT	1	30041102
EasyPlus drying tube	1	30044701
Valve	1	30042860
Bottle head	1	30060023
Insert set titration head GT	1	30057635
Burette cover set	1	30057633

Accessories

Product	Order Number
EasyPlus burette 10 mL	30043900
EasyPlus burette 50 mL	30043902
EM43-BNC Sensor Redox Ipol	30043105
Bottle adaptor China	30060026
Bottle adaptor India	30060024
Bottle adapter Merck	23774
Bottle adaptor Fischer	23787
EasyPump	30064218
Siphon tips (5 pcs.)	23240
Glass beakers (80 mL), 20 pcs.	101446
Plastic beaker (80 mL), 120 pcs.	51109388
EasyDirect instrument license for 3 titrators per PC	30065449
USB-P25 strip printer	11124301

Technical specifications

mV sensor input	Measuring range	± 2000 mV
	Resolution	± 0.1 mV
Polarized sensor input	Current source redox ipol	0 – 5 µA
	Resolution	± 0.1 µA
	Measuring range (voltage input)	0 – 2000 mV
Power supply	Input voltage	24 VDC / 1.25 A
	Line input voltage	100 – 240 V ±10 %
	Input frequency	50/60 Hz
Environment	Ambient temperature	5 °C – 40 °C
	Atmospheric humidity	Max. 80% at 31 °C
Titration dimensions	Width x depth x height / weight	170 x 220 x 350 mm / 1.9 kg (Plus KF stand 0.8 kg)
Control and display	Touchscreen	4.3 inch, 480 x 272 pixel RGB
Communication	USB-A	USB full / low speed
	USB-B	USB full / low speed
Application*	Repeatability	0.5% RDS
	Titration type	Direct, back compensated, blank compensated

*Depends on the environment and handling

Support or Service? No problem



Application know-how

Would you like application know-how to titrate a specific sample?



Expertise & help

Would you like best practice titration tips or require technical support expertise?



Maintenance & verification

Do you want learn more about preventive maintenance or need to carry out a performance verification of your titrator?

Simply register to the EasyPlus™ Titration Webpage and quickly find comprehensive help for your daily tasks.

► www.mt.com/easyplustitration

METTLER TOLEDO Group

Analytical

Local contact: www.mt.com/contact

Subject to technical changes.

© 06/2018 METTLER TOLEDO

30073712A

Marketing Titration / MarCom Analytical

www.mt.com/easyplustitration

For more information