

Material Number: 85701.601

Material description: Ascorbic acid test strips, Dosatest

Date specification set: 31st Mar 2021

(All batches of this product after this date are tested to this specification. In the unlikely event that the specification is changed a new certificate will be issued)

Shelf life (from production)

2.5 Years

Scope of delivery

100 test strips in a tube, instruction leaflet

Parameter

Ascorbic acid

Gradation

0-50-100-200-300-500-700-1000-2000 mg/L Vitamin C

Method

Phosphomolybdenum blue

Color change

Yellow → green-blue

Signature

This document has been produced electronically and is valid without a signature.

Signed on behalf of VWR:

Ben Petri

Quality Manager