

Technical Data Sheet



VWR® 10mL Open End Serological Pipettes

Description

10mL in 1/10mL

Purpose

Use for liquid transfer


Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI

Color: Clear

Filter: PP (Polypropylene)

Color: White



Unit: mm

	Dimensions	Tolerance
A	ø10.0	±0.1
B	ø8.2	±0.1
C	ø7.95	±0.35
D	308.5	±2.0

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: orange
- Graduation (mL): in 1/10mL
- Negative graduations up to -3 ml (total volume = 13 ml)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 10mL Open End Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5861	76292-294	10mL, Orange, Individually Wrapped, Paper/Plastic	Yes	1	/	400
612-5862	76292-296	10mL, Orange, Bulk packed	Yes	/	25	400