

Material Number: 87042.180
Material description: Ion chromatography cation multi-element standard solution VI
Grade: ARISTAR®
Date specification set: 25th Sep 2019

(All batches of this product after this date are tested to this specification. In the unlikely event that the specification is changed a new certificate will be issued)

Characteristics	Specifications	Traceability
NH ₄ ⁺	9.98 – 10.02 mg/l	NIST SRM No 194a
K ⁺	49.9 – 50.1 mg/l	NIST SRM No 3141a
Na ⁺	99.8 – 100.2 mg/l	NIST SRM No 3152a
Ca ²⁺	99.8 – 100.2 mg/l	NIST SRM No 3109a
Mg ²⁺	99.8 – 100.2 mg/l	NIST SRM No 3131a

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International