

Certificate of Conformity

VWR Catalogue Number	See List
Description	Disposable cuvettes with 10 mm Path Length
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2025-07-11

Catalogue Number	Type	Range [nm]	Volume [mL]	Window Material
634-0678	Semi-micro	300 – 900	1.5 – 3	PMMA
634-0677	Macro	300 – 900	2.5 – 4.5	PMMA
634-0676	Semi-micro	340 – 900	1.5 – 3	PS
634-0675	Macro	340 – 900	2.5 – 4.5	PS

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material Polystyrene (PS) or Polymethylmethacrylate (PMMA). The external dimension of disposable cuvettes macro, semi-micro and micro are manufactured according to DIN 58963 part 3.

Polystyrene Cuvettes: Standard deviation at 360nm $\leq \pm 0.005$ extinction units

Polymethylmethacrylate Cuvettes: Standard deviation at 320nm $\leq \pm 0.004$ extinction units

Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	A small quantity of a substance from animal origin (tallow derived substance) was introduced during the manufacturing process.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethylhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Should be stored using the original packaging, at temperature between 10°C and 50°C. Relative Humidity between 30% and 65%.
Shelf Life	5 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	12-2021
2	2.0	Modified Quality system compliance and disclaimer	07-2025