

11/08/2022

Freudenberg Medical Europe
2-8 Liebigstrasse
67661 KAISERSLAUTERN
GERMANY

PO# 4500030335

QTY. 2 EA

≈ 50 St.

PRODUCT CERTIFICATION

FREUDENBERG P/N DESCRIPTION FREUDENBERG L/N

70-400-07

0000231470

**¾" SILICONE GASKET
DOW CORNING, SILASTIC Q7-4780 MATERIAL**

MATERIAL FM P/N FM L/N

Silicone	30-011-01	0000211019
Elastomer	30-012-01	0000211020

Freudenberg Medical certifies that our silicone gaskets are manufactured following Good Manufacturing Practices (GMP) in accordance with the new Quality System Regulation (QSR) in a FDA registered medical device facility. Freudenberg Medical has been certified to meet ISO 13485 quality system. Material manufacturer states that these materials meet or exceed USP Class VI biocompatibility requirements and FDA Regulations CFR 177.2600 for contact with food and drugs. These platinum-cured silicone products have full traceability to raw materials and contain no peroxide catalysts.

Sincerely,

Benjamin Echevarria / Shipping Clerk

Freudenberg Medical	
Material Number	Batch Number
45633007	0000231470

(7040007 REV. B)

11/08/2022

Freudenberg Medical Europe
2-8 Liebigstrasse
67661 KAISERSLAUTERN
GERMANY

PO# 4500030335

QTY. 1 EA

25 Hz

PRODUCT CERTIFICATION

FREUDENBERG P/N DESCRIPTION FREUDENBERG L/N

70-400-07


0000234331

**¾" SILICONE GASKET
DOW CORNING, SILASTIC Q7-4780 MATERIAL**

<u>MATERIAL</u>	<u>FM P/N</u>	<u>FM L/N</u>
Silicone	30-011-01	0000211019
Elastomer	30-012-01	0000211020

Freudenberg Medical certifies that our silicone gaskets are manufactured following Good Manufacturing Practices (GMP) in accordance with the new Quality System Regulation (QSR) in a FDA registered medical device facility. Freudenberg Medical has been certified to meet ISO 13485 quality system. Material manufacturer states that these materials meet or exceed USP Class VI biocompatibility requirements and FDA Regulations CFR 177.2600 for contact with food and drugs. These platinum-cured silicone products have full traceability to raw materials and contain no peroxide catalysts.

Sincerely,

Benjamin Echevarria / Shipping Clerk

Freudenberg Medical	
Material Number	Batch Number
45639007	0000234331

(7040007 REV. B)