



Specification

1.09877.0010 OGYE selective supplement

Specification		
Tested in OGYE-agar base art. no. 1.05978.		
Recovery on test medium (Candida albicans ATCC 10231 (WDCM 00054))	≥ 50	%
Recovery on test medium (Saccharomyces cerevisiae ATCC 9763 (WDCM 00058))	≥ 50	%
Recovery on test medium (Aspergillus brasiliensis (formerly A. niger) ATCC 16404 (WDCM 00053))	≥ 50	%
Recovery on test medium (Penicillium aurantiogriseum ATCC 16025)	≥ 50	%
Growth (Bacillus subtilis ATCC 6633 (WDCM 00003))	none	
Growth (Escherichia coli ATCC 25922 (WDCM 00013))	none	

Incubation: 3 up to 5 days; 28 °C; aerobic.

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.