

Corning® Centrifuge Tubes

Designed for Optimum Performance

Corning offers a full line of disposable, premium quality conical centrifuge tubes manufactured from high quality, ultra-clear resins. These tubes meet the quality and performance standards demanded by researchers.

► High Speed 15 mL and 50 mL tubes

- Rated up to 17,000 RCF
- 95 kPa (14 psi) pressure-tested (plug seal cap)
- RNase-/DNase-free
- Nonpyrogenic
- Sterile

► Large Volume 250 mL and 500 mL tubes

- Rated up to 6,000 RCF
- Best choice for cell harvesting and purification
- Nonpyrogenic
- Sterile



Ordering Information

Corning® 15 mL Centrifuge Tubes



15 mL and 50 mL centrifuge tubes with Corning CentriStar cap

VWR Cat. No.	Corning Cat. No.	Description	Material	Cap Style	Max. RCF	Qty/Pk	Qty/Cs
21008-667	430053	Conical bottom	PET	Plug seal	3,600	50/sleeve	500
21008-656	430055	Conical bottom	PET	Plug seal	3,600	50/rack	500
21008-678	430052	Conical bottom	PP	Plug seal	12,000	50/rack	500
21008-673	430766	Conical bottom	PP	Plug seal	12,000	50/sleeve	500
21008-670	430790	Conical bottom	PP	CentriStar™	12,500	50/rack	500
21008-671	430791	Conical bottom	PP	CentriStar	12,500	50/sleeve	500

Corning 50 mL Centrifuge Tubes



250 and 500 mL centrifuge tubes with support cushions

VWR Cat. No.	Corning Cat. No.	Description	Material	Cap Style	Max. RCF	Qty/Pk	Qty/Cs
21008-775	430897	Self-standing bottom	PP	Plug seal	3,000	25/sleeve	500
21008-777	430921	Self-standing bottom	PP	CentriStar	3,000	25/sleeve	500
21008-758	430304	Conical bottom	PET	Plug seal	3,600	25/rack	500
21008-714	430290	Conical bottom	PP	Plug seal	15,500	25/rack	500
21008-725	430291	Conical bottom	PP	Plug seal	15,500	25/sleeve	500
21008-690	430828	Conical bottom	PP	CentriStar	17,000	25/rack	500
21008-736	430829	Conical bottom	PP	CentriStar	17,000	25/sleeve	500
89089-818	4558	Conical bottom	PP	CentriStar	17,000	25/rack	300

250 mL and 500 mL Premium Quality Centrifuge Tubes and Support Cushions

VWR Cat. No.	Corning Cat. No.	Description	Material	Cap Style	Max. RCF	Qty/Pk	Qty/Cs
21008-771	430776	250 mL Tube	PP	Plug seal	6,000	6	102
21008-780	430236	250 mL Support cushion*	PEI	N/A	N/A	N/A	6
89091-000	431123	500 mL Tube	PP	Plug seal	6,000	6	36
29443-019	431124	500 mL Support cushion*	PEI	N/A	N/A	N/A	6

*Support cushions must be used with the 250 mL and 500 mL centrifuge tubes unless the rotor has appropriately shaped V-bottom holders. PP = Polypropylene, PET = Polyethylene Terephthalate, RCF = Relative Centrifugal Force (x g), PEI = Polyetherimide

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your VWR representative.


800 932 5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.