

Bucket, stainless steel



The bucket is made entirely of high-gloss polished stainless steel. The smooth surface of the stainless steel bucket makes it easy to clean. The rounded, movable handle makes it comfortable to move the bucket. The laboratory bucket is suitable for moving liquids and bulk materials properly, especially in the chemical, pharmaceutical or cosmetic industries, and in the food and animal feed industries.

- Stainless steel (1.4301)
- High-gloss polished

Bucket, stainless steel					
Capacity l	Dia. mm	Height mm	Item no.		
10	275	275	2309-0010	N4	
15	330	305	2309-0015	N4	

Packaging bucket



Packaging bucket for transport, packaging and storing of various materials. Suitable for the food or building material industries and non-aggressive chemicals in environmental protection, hospitals and the household. With strong handle. Quality closure on the bucket for product safety.

- PP, white
- Air-tight closure
- Complete with lid
- With advertising printing in large quantities on request

Bulk prices

Packaging bucket					
Capacity l	Dia. top mm	Height with cap mm	Item no.		
1	132	131	2327-0001	N3	
2	175	132	2327-0020	N3	
3	225	123	2327-0030	N3	
5	225	194	2327-0050	N3	
10	268	263	2327-0010	N3	