

Werksbescheinigung / Declaration of compliance with the order

gemäß / in accordance with DIN EN 10204 - 2.1

Mischgefäß / Tube:

Modell / Model:

ST-20/50..., DT-20/50..., BMT-20/50...
TC-20-M, TC-50-M
DIS 50-T-M, DIS 50 S-M, DIS 300-S, DIS 300-S-M

Produktberührende Materialien / Parts in contact with medium.

Deckel / Cover:	Kunststoff / Plastic: PP
Deckel / Cover: (DIS-300-S)	Kunststoff / Plastic: PE
Deckel / Cover: (DIS-50-S-M, DIS-50-T-M)	Kunststoff / Plastic: PP / TPE
Membran / Membrane:	Kunststoff / Plastic: PP
Probengefäß / Testing vessel:	Kunststoff / Plastic: PP
Rührorgan / Stirring stick:	Kunststoff / Plastic: PEEK
Rührorgan / Stirring stick: (DIS-300-S, DIS-300-S-M)	Kunststoff / Plastic: PEEK Compound
Dispersiereinheit / Dispersing unit:	Kunststoff / Plastic: PEEK Compound
Dissolver / Dissolver: (DIS-50-S-M, DIS-300-S, DIS-300-S-M)	Edelstahl / Stainless steel: 1.4404
Dissolver / Dissolver: (DIS-50-T-M)	Titan / Titanium: 3.7165
Glaskugel / Glass ball: (BMT-20/50)	Glas / Glass: Kalknatronglas / Soda lime glass:
Metallkugel / Metal ball: (BMT-20/50)	Edelstahl / Stainless steel: (1.4301) 1.4401

Staufen, 25.11.2021

20000012565 D

Stefan Meidroth
Director
Quality Assurance

Hans-Urs Eckerle
Team Leader
Regulatory Compliance

Dieses Dokument wurde elektronisch erstellt und ist ohne Unterschrift gültig. (DIN EN 10204, Abschnitt 5) /
This document was created electronically and is valid without signature. (DIN EN 10204, Chapter 5)



designed
to work perfectly

Brazil	India	Malaysia	China	USA	Germany
IKA® do Brasil Rua Matão, 122 Campinas-SP CEP 13030-300, Brazil South and Central America	IKA® India Private Limited 814/475, Survey No. 129/1 Mysore Road, Kengeri, 560060 Bangalore, Karnataka India	IKA® Works (Asia) Sdn Bhd Lot PT6445, Jalan Industri 3/4 Rawang Integrated Industrial Park 48000, Rawang, Selangor, Malaysia	IKA® Works Guangzhou 173 – 175 Friendship Road Guangzhou Economic and Technological Development District 510730 Guangzhou, China	IKA® Works, Inc. 2635 Northchase Pkwy SE Wilmington, NC 28405-7419 USA	IKA®-Werke GmbH & Co.KG Janke & Kunkel-Str. 10 79219 Staufen Germany
Tel. +55 19 3772-9600 Fax +55 19 3772-9601 info@ika.net.br	Tel. +91 80 26253900 Fax +91 80 26253901 info@ika.in	Tel. +60 3 60995679 Fax +60 3 6092 0240 sales@ika.my	Tel. +86 20 8222-6771 Fax +86 20 8222-6776 info@ika.cn	Tel. +1 910 452-7059 Fax +1 910 452-7693 sales@ika.net	Tel. +49 7633 831-0 Fax +49 7633 831-98 sales@ika.de