



Certificate of Conformity

VWR TEST SIEVES, 76x22 mm

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	Ø (mm)	Height (mm)	Mesh Size (mm)	Pk
510-3611	76	22	0.020	1
510-3612	76	22	0.025	1
510-3613	76	22	0.032	1
510-3614	76	22	0.036	1
510-3615	76	22	0.038	1
510-3616	76	22	0.040	1
510-3617	76	22	0.045	1
510-3618	76	22	0.050	1
510-3619	76	22	0.053	1
510-3620	76	22	0.056	1
510-3621	76	22	0.063	1
510-3622	76	22	0.071	1
510-3623	76	22	0.075	1
510-3624	76	22	0.080	1
510-3625	76	22	0.090	1
510-3626	76	22	0.100	1
510-3627	76	22	0.106	1
510-3628	76	22	0.112	1
510-3629	76	22	0.125	1
510-3630	76	22	0.140	1
510-3631	76	22	0.150	1
510-3632	76	22	0.160	1
510-3633	76	22	0.180	1
510-3634	76	22	0.200	1



Certificate of Conformity

Catalog No.	Ø (mm)	Height (mm)	Mesh Size (mm)	Pk
510-3635	76	22	0.212	1
510-3636	76	22	0.224	1
510-3637	76	22	0.250	1
510-3638	76	22	0.280	1
510-3639	76	22	0.300	1
510-3640	76	22	0.315	1
510-3641	76	22	0.400	1
510-3642	76	22	0.425	1
510-3643	76	22	0.450	1
510-3644	76	22	0.500	1
510-3645	76	22	0.560	1
510-3646	76	22	0.600	1
510-3647	76	22	0.630	1
510-3648	76	22	0.710	1
510-3649	76	22	0.800	1
510-3650	76	22	1.00	1
510-3651	76	22	1.12	1
510-3652	76	22	1.18	1
510-3653	76	22	1.25	1
510-3654	76	22	1.40	1
510-3655	76	22	1.60	1
510-3656	76	22	1.70	1
510-3657	76	22	1.80	1
510-3658	76	22	2.00	1
510-3659	76	22	2.24	1
510-3660	76	22	2.36	1
510-3661	76	22	2.50	1
510-3814	76	22	0.900	1
510-3815	76	22	0.850	1
510-3816	76	22	0.315	1



Certificate of Conformity

Material

Frame and Tension ring = AISI 304 (1.4301)

Test Sieve Cloth = AISI 316 (1.4401) for Mesh sizes 0.020 – 2.50 mm

Test Sieve Cloth = AISI 316 (1.4301) for Mesh sizes 2.80 – 125 mm

Compliance

- In compliance to DIN ISO 3310-1, ISO 3310-1, ASTM E11, BS ISO 3310-1
- Maintenance of Test Sieve is in compliance to ISO 2591-1, DIN 66165

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2017**