



Certificate of Conformity

VWR SURFACE SAMPLING SPONGES

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

VWR Cat. No.	Description	Dimensions & volume	Closure	Volume of solution
300-0229	Sponge dosed with buffer	177.80 mm x 304.80 mm	2 round wires	6.5 ml
300-0230	Surface sampling sponge with neutraliser	139.70 mm x 304.80 mm	Open top	10 ml
300-0231	Surface sampling sponge with disinfection	139.70 mm x 304.80 mm	Open top	10 ml
300-0232	Surface sampling sponge with handle	139.70 mm x 304.80 mm	2 round wire closure system	10 ml
300-0240	Surface sampling sponge with handle	139.70 mm x 304.80 mm	2 round wire closure system	10 ml

Closure Material: Steel wires covered with FDA approved pressure sensitive tape with a backing of polypropylene or PVC.

Thickness of Material : 3 ml

Compliance: Comply to the applicable regulations of the Food and Drug Administration (21 CFR 177.1520 (c)3.1a, (c)3.2a and 178.2010) for use in contact with food.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2015**