

VWR® Slide Zipper bags



Description: VWR® Slider Zipper bag, LDPE, transparent, with blue slider

Purpose: Packing goods

Materials: LDPE 260GG+FD21SH

Features:

- Uses: under -10 ~65 °C
- Sterilization: N/A
- BSE/TSE: The above product contains no ingredients of animal origin and no material derived
 - from, or exposed to animals affected by or under quarantine for Transmitting, Animal
 - Spongiform Encephalopathy (TSE)/Bovine Spongiform Encephalopathy (BSE).
- Products are manufactured under the ISO 9001:2015 standard.
- Shelf Life: N/A.
-

VWR® Slider Zipper Bag

EU Part No.	Product Description		
VWRI129-1202	VWR® Slider Zipper bag, LDPE, transparent, 250mmx350mm+92Bottom Gussetx50um , with blue slider		
VWRI129-1203	VWR® Slider Zipper bag, LDPE, transparent, 300mmx400mm+92Bottom Gussetx50um, with blue slider		

EU Part No.	Width(mm)	Length(mm)	Thickness(um)	Printing	Printing size	Side strength	Tensile strength		Elongation
	Tolerance(+3%/-2%)	Tolerance(+3%/-2%)	Tolerance(±5%)				MD	TD	
VWRI129-1202	250	350+92BG	50	NO	N/A	10N/15mm	21	21	550%
VWRI129-1203	300	400+92BG	50	NO	N/A	10N/15mm	21	21	550%