

Sodium Hydroxide
Baker Analyzed
Pellets



Material No: 0402.5000
Effective from: 2019/06/05

Product Specification

Test	Specification
Assay	$\geq 98 \%$
Carbonate (as Na_2CO_3)	$\leq 1 \%$
Potassium (K)	$\leq 0.05 \%$
Silica (SiO_2)	$\leq 0.001 \%$
Trace Impurities (in ppm)	
Aluminium (Al)	≤ 5
Calcium (Ca)	≤ 5
Chloride (Cl)	≤ 60
Heavy Metals (as Pb)	≤ 5
Iron (Fe)	≤ 5
Nickel (Ni)	≤ 5
Phosphate (PO_4)	≤ 5
Sulphate (SO_4)	$\leq 30 \text{ ppm}$
Total Nitrogen (N)	≤ 3

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000