



Laboratory glassware,
apparatus and
precision instruments



**Glaswarenfabrik
Karl Hecht
GmbH & Co KG**

Stettener Str. 22 - 24
97647 Sondheim v. d. Rhön
Germany

Tel. +49 9779 808-0
Fax +49 9779 808-88

Email: info@assistent.eu
Internet: <http://www.assistent.eu>
VAT No. DE 132 198 535

Data sheet

Date of issue: 14.11.2016

Revised: 08.08.2019

Product description

Test tubes, AR-Glas® (soda glass), heavy wall, round bottom, without rim

	Nominal size	Wall thickness	also known as
4 2775 003	50 x 6 mm	0.65 mm	
4 2775 006	27 x 7 mm	0.7 mm	
4 2775 009	35 x 6 mm	0.65 mm	Durham tubes
4 2775 012	40 x 8 mm	0.7 mm	
4 2775 015	80 x 10 mm	0.9 mm	
4 2775 018	70 x 10 mm	0.9 mm	
4 2775 021	100 x 10 mm	0.9 mm	
4 2775 024	75 x 12 mm	0.7 mm	Kahn tubes
4 2775 027	100 x 12 mm	0.7 mm	tubes for blood typing
4 2775 030	100 x 14 mm	0.9 mm	
4 2775 033	130 x 14 mm	0.9 mm	
4 2775 036	75 x 15 mm	0.9 mm	
4 2775 039	85 x 15 mm	0.9 mm	
4 2775 042	135 x 15 mm	0.9 mm	
4 2775 045	100 x 16 mm	0.8 mm	
4 2775 051	160 x 16 mm	0.8 mm	
4 2775 054	100 x 18 mm	1.0 mm	
4 2775 060	180 x 18 mm	1.0 mm	
4 2775 251	160 x 16 mm	0.8 mm	

Packing: box of 100 - except No. 4 2775 024 and 4 2775 251: box of 200

Manufacturing process

The test tubes are of automatic production. They are made of soda glass tubing.

The cutting of the tubing as well as the shaping of the bottom and fusing of the rim are automatically done. At the end of the production process, the test tubes are thermally treated.

Photos

