

PRODUCT SPECIFICATION /CERTIFICATE OF CONFORMITY

www.dwkLtd.com
enquiries.uk@dwk.com

Specification	Product Code	JPB2000P
	Description	Azlon 2000ML Jug with handles, blue printed graduations.
	Material	Polypropylene
	Capacity	2000
	Colour	Natural
	Autoclavable	Yes
	Pack quantity	1

Dimension	Height (mm)	214
	Width (mm)	150
	Rim diameter (mm)	150
	Graduation (ml)	500

The above product is compliant to the following standards only if marked 'Yes'

Compliance	Food safe (EU/FDA)	EU 10/2011 / FDA CFR 21 177
	BSE / TSE free	YES
	REACH Compliant	YES
	LATEX	YES

DWK Life Sciences Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of DWK Life Sciences Limited by:

P Burton
Quality Technician

12 Jan 2021