

Certificate of Conformity

VWR Catalogue Number	See List
Description	Test tubes with red cap, sterile
Country of Origin	Manufactured in Italy
Date of Issue (yyyy-mm-dd)	2021-10-20

Catalogue Number	Capacity [ml]	Length [mm]	External Diameter [mm]	Material
216-1166	4	75	13	Polypropylene
212-0726	5	75	13	Polystyrene
216-1047	10	100	16	Polypropylene
216-1045	10	100	16	Polystyrene

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

Material

See above table.

The products are considered safe to use in food, pharmaceutical and cosmetic contact applications.

Sterilization

Products are sterilized using Beta rays and meets SAL 10⁻⁶ with a specified dosage range of 17 – 20 kG.

ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
GMO Statement	Genetically modified materials (GMO) have not been intentionally incorporated in the manufacturing of the product.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The products are not made from natural rubber latex.
BPA Statement	Bisphenol A (BPA) has not been intentionally incorporated in the manufacturing of the products.
DEHP Statement	DEHP has not been intentionally incorporated in the manufacturing of the products.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store in a dry place, away from sources of heat and direct sunlight.
Shelf Life	60 months.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature