



World Class Quality
ISO 9001:2015

PRODUCT SPECIFICATIONS

PRODUCT #331 - 90x15mm

LID: VENTS, NO RIBS
BASE: 3 SECTION WITH "ISO PLATER" FROST MARK

Raw Material	• 100% virgin, high molecular, crystal clear polystyrene • Produced by an ISO 9001 registered company		
Weight of one Complete Assembly (g)	12.60 $\pm$ 0.30		
Height with stacking rings (mm)	13.58 $\pm$ 0.20		
Divider (mm)	HIGHER SIDE: 10.55 $\pm$ 0.20 LOWER SIDE: 8.65 $\pm$ 0.20		
Stackability	Beveled stacking ring for stability "S e m i - S t a c k i n g"		
Diameter (mm)		OUTSIDE	INSIDE
	BASE-TOP	87.50 $\pm$ 0.2	85.20 $\pm$ 0.2
	BASE-BOTTOM	85.00 $\pm$ 0.2	84.10 $\pm$ 0.2
	LID-TOP	89.50 $\pm$ 0.2	88.50 $\pm$ 0.2
	LID-BOTTOM	92.75 $\pm$ 0.2	89.10 $\pm$ 0.2
Ventilation Lugs (mm)	QTY.: 4 HEIGHT: 0.25 $\pm$ 0.025		
Flatness (mm)	BASE: +0.25 to -0.25 LID: +0.20 to -0.50		
Transparency	Excellent		
Mode of Production and Sterility	CLEAN ROOM "ISO Class 5" and Phoenix own controlled environment cleanliness limits		
Packaging	25 dishes per sleeve (HDPE 2.5 mil) 24 sleeves per carton with 2 inner liners		
Production Number (LOT# & EMNS)	LOT# MDDPPN (Month (A,B,C...), Day, Phoenix Product Number) EMNS (Employee Number, Shift)		
Shelf Life	84 MONTHS EXP.DATE YYYY-MM (Year-Month)		
Certificate of Sterility	YES Sterility testing in accordance with current USP procedure, conducted by an independent Microbiological Laboratory.		

Rev. # 9

PHOENIX BIOMEDICAL PRODUCTS S.L

Ctra. El Estrecho-Lobosillo Km.2, Parque Tecnológico de Fuente Álamo
Avda. Emilio Pérez Piñero, Parcela 3.7
30320 Fuente Álamo, MURCIA (SPAIN)

Tel: (34) 968 197610 • Fax: (34) 968 197602

E-mail: stardish@phoenix-biomed.com

www.phoenix-biomed.com