

BAXTER CLEVELAND PRODUCT CERTIFICATION

**911 North Davis
Cleveland, MS 38732**

(2F7121 / 2F7122 / 2F7123 / 2F7124 / 2F7125)

CODE <u>2F7124</u>	BATCH <u>0144216</u>	MFG. DATE <u>3/14/20</u>	STERILIZATION DATE <u>3/14/20</u>
PLANT OF MANUFACTURE CLEVELAND, MS		EXPIRY DATE <u>March 2023</u>	

DESCRIPTION

0.9% Sodium Chloride, USP

TEST	LIMITS	SAMPLE RESULTS	DATE PASSED
BACTERIAL ENDOTOXIN (PYROGEN)	CONTROL LIMIT (Eu/mL)		
13-01-10-313	<0.5 EU/ml = Pass	<u><0.01</u> EU/ml	<u>3/16/20</u>
STERILITY PARAMETRIC RELEASE	PASS	PASSED	<u>3/14/20</u>
PARTICLE ANALYSIS 20-03-07-006	NMT 25 part/mL $\geq 10 \mu\text{m}$	<u>0</u> part/mL	<u>3/16/20</u>
	NMT 3 part/mL $\geq 25 \mu\text{m}$	<u>0</u> part/mL	

PARAMETRIC RELEASE IN ACCORDANCE WITH U.S. FDA REQUIREMENTS

CHEMICAL RESULTS 3-15-19-156 ISSUE DATE 11/13/18

REFERENCE	TEST	LIMITS	RESULT
11-21-16-003	pH @ 25°C	4.5-7.0	<u>5.5</u>
11-07-19-001 or 11-16-03-005 or 11-16-03-016	NaCl Assay	8.55-9.45 g/l	<u>8.96</u> g/l
11-20-19-006	Sodium ID - USP	Positive	Positive

THIS PRODUCT WAS MANUFACTURED ACCORDING TO CURRENT GMP, FULLFILLS THE REQUIRMENTS OF THE MASTER PRODUCTION DOCUMENT AND TESTED TO USP COMPLIANCE.

DATE BATCH RELEASED 3/14/20PREPARED BY/DATE Brooklynson 3/14/20QUALITY MANAGEMENT SUPERVISOR [Signature] 3-19-20

MDM

1. Code, Batch, Mfg. Date, Sterilization Date, and Expiry Date - Enter the information from Production Report of batch required.

2. Test:

Bacterial Endotoxin – Enter the bacterial endotoxin test results. Enter the date passed (date the bacterial endotoxin testing was complete) from the bacterial endotoxin test report for the batch.

Sterility – Enter the date product was sterilized.

Particle Analysis – Enter the sample results and the date passed from the particle analysis record for the batch.

3. Chemical Results - Enter issue date of the applicable 3-15 specification. All results should be verified against the limits in the applicable 3-15 specification. All test results are retrieved from the EOB chemical testing. Fill in test results.

Note: Multiple test results may be averaged by Release Clerk and the average entered on this certification.
(Example - multitank and stability batches).

Note: NaCl Assay- Circle the 11 spec used for testing.

4. Date Batch Released - Fill in date batch was released.

5. Prepared By/Date - Dated signature of person entering information.

Note: Signature of Release Clerk indicates the results have been verified against the results listed on the chemical final analysis and the applicable 3-15 chemical testing specification, the results listed on the Particle Analysis Report and results listed on the Bacterial Endotoxin report.

6. N/A sections not used.

7. Quality Management Supervisor - Signature of Quality Management Supervisor or delegate verifying the Batch #, Code, Manufacturer Date, Expiration Date, and Batch Release Date documented matches the information documented on the Production Record. Individual will and verify that all test results listed are within specification and/or SOP limits.

DESCRIPTION OF CHANGE		
Update MDM to include specific information for the approver to verify.		
REASON FOR CHANGE		
Comply with CL-03-01-039.		
CHANGE HISTORY		
CHANGE #	ISSUE DATE	DESCRIPTION OF CHANGE
19-0877	04/12/19	13-1-V is obsolete and adding 11-21-16-005 test to Particle Analysis. Eu/mL is being added to sample results for Bacterial Endotoxin and part/mL is being added to Particle Analysis sample results. 11-20-19-007 is being deleted because it's not used for these codes. Add change history to form.
N/A	N/A	N/A
N/A	N/A	N/A
APPROVALS ON FILE		
ROLE	NAME	
SME/Release Clerk	Beverly Robinson	
Approver	Jayme Williams	
Concurrer	Ellen Soliz	