

Certificate of Conformity

VWR Catalogue Number	See List
Description	Microscope slides VWR®, Clear, white Glass, ISO 8037-1 certificated
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-09-09

Catalogue Number	Designation	Thickness [mm]	Length x Width [mm]
631-1553	Slide with bevelled edges and frosted end	1.0	76 x 26
631-1552	Slide with bevelled edges and plain end	1.0	76 x 26
631-1551	Slide with cut edges and frosted end	1.0	76 x 26
631-1550	Slide with cut edges and plain end	1.0	76 x 26

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [ISO 9001:2015](#).

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Float Glass.
Refractive index	1.53 +/-0.02 of refractive index of the light at the green mercury e-line (wavelength=547.07nm) ISO 8037-2.
Wettability	The central parts of the two surfaces of the slides should be wetted equally, and slightly dehydration of the edges is allowed. (ISO 8037-2)
Sterilization	Not applicable.

ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Room temperature, dry.
Shelf Life	No shelf life.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.


VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.


Annex

Technical Drawing


631-1550




631-1551



631-1552



631-1553



Version History

Sl.no	Version number	Description	Month & Year
1	1.1	Created General CoQ using new layout	09-2022
2	2.0	Modified Quality system compliance and disclaimer	09-2025