

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Spinwedge® Magnetic Stir Bars
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-04-20

VWR Catalogue Number	Diameter	Length
76001-896	12.7 mm (1/2")	9.5 mm (3/8")
58949-232	25 mm (1")	9.5 mm (3/8")
58949-254	44 mm (1 3/4")	12 mm (1/2")
58949-265	50 mm (2")	12 mm (1/2")

Quality System Compliance

Products are manufactured under the [ISO 9001:2015 & ISO/IEC 17025:2017](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

- Tolerance +/- .020 "
- Autoclavable @121° C; for 20 minutes

Product Specifications

Material	Dupont White Teflon #807N-X; Magnet, ALNICO 5; Materials Meet FDA Regulations 21CFR177.1550, 21CFR177.1520 Teflon used to totally encapsulate the magnetic bar meets US Pharmacopeia's Class VI Based on the information provided to us this product does not contain or come in contact with lead, latex, Mercury, Formaldehyde or other contaminants.
Sterilization	No test performed.
ATP Assay	No test performed.
DNase & RNase Free	No test performed.

BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	No test performed.
Non-Pyrogenic Statement	No test performed.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Store at ambient warehouse temperatures, low humidity, protected from direct sunlight.
Shelf-Life	Product has no expiry date.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.