

Material No.: 16004-382

Material Description: VWR Microscope Slide, Unimark, 75x25x1mm, White, 72/box, 20 boxes/cs

Lot No.: 54818

Compliance and Procedure Result

Product Specifications

- VWR VistaVision Unimark Microscope Slide, 72/Box, 20 Boxes/cs.
- Length x Width: 75.2 x 25.2 mm
- Thickness: 1.00 mm
- Country of origin: Germany
- Soda lime glass
- Pre-cleaned, edges ground at 90°, white marking area of approx. 20mm on one side
- Ready for use

Accuracy

- Tolerance of dimensions: +/- 0.2 mm
- Tolerance of thickness: +/- 0.05 mm

Compliance

Standard: DIN ISO 8037-1

Signed on behalf of VWR By: Janhvi Nanda

Quality Employee

7/29/2025