

Certificate of Conformity

Catalogue Number	49002-01061
Description	CONDUCTIVE TIPS 200 µL FILTER TRAY CS2304
Lot no.	1000121/1-1
Country of Origin	Manufactured in Germany
Date of Issue (dd/mm/yyyy)	27/10/2021
Production date	2021/10
Expiry date	2026/10

Quality System Compliance

Products are manufactured in [ISO 13485:2016](#) and [ISO 9001:2015](#) certified plant. Products are inspected and controlled in accordance with the current applicable product specification.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Tips: Polypropylene Tray: Polycarbonate or For tips H, R or C-tray: Polypropylene
Manufacturing Environment	Robotic tips are manufactured under clean room conditions according to DIN EN ISO 14644-1:2016-06.
Sterilization	Not applicable
ATP Assay	Not applicable
PCR inhibitors Free	Product yield equivalent to control < 10 amplifiable targets
Human DNA Free	Tested with detection level < 1 pg
DNase & RNase Free	Tested with detection level < 1*10 ⁻⁶ KunitzU (DNase) and < 1*10 ⁻⁹ KunitzU (RNase)
Nucleic Acid Free	Not applicable
BSE/TSE	No use of any raw materials produced from or substances derived from animal origin. The manufacturing process of the product does not intentionally use any ingredient of animal origin and no material

derived from or exposed to animals affected by.

Cytotoxicity	Not applicable
Non-Pyrogenic Statement	Not applicable
Latex Statement	Not applicable
BPA Statement	Not applicable
DEHP Statement	Not applicable
RoHS	The raw materials used are compliant with the currently valid directive 2011/65/EU and 2002/95/EG (RoHS).
REACH Statement	The product does not contain any "SVHC" ingredients, in accordance with the currently valid regulation (EC) No 1907/2006 on REACH.
Material Statement	Not applicable
Storage Conditions	Keep dry and away from sunlight
Shelf Life	5 years

Disclaimer: Ritter states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.