



Specification

1.05448.0500 Wort agar for microbiology

Specification

Appearance (clearness)	clear
Appearance (colour)	brownish
Solidification behaviour (2 hrs., 45 °C)	liquid
pH-value (25 °C)	4.8 - 5.2

Specification

Recovery on test medium (Candida albicans ATCC 10231 (WDCM 00054))	≥ 70	%
Recovery on test medium (Saccharomyces cerevisiae ATCC 9763 (WDCM 00058))	≥ 70	%
Recovery on test medium (Saccharomyces cerevisiae ATCC 9080)	≥ 70	%
Recovery on test medium (Rhodotorula mucilaginosa DSM 70403)	≥ 70	%
Growth (Geotrichum candidum DSM 1240)	good to very good	
Growth (Penicillium commune ATCC 10428)	fair to very good	
Growth (Aspergillus brasiliensis (formerly A. niger) ATCC 16404 (WDCM 00053))	good to very good	
Growth (Trichophyton ajelloi ATCC 28454)	fair to good	

Incubation: 7 days; 28 °C; aerobic (yeasts 3 days).

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.