


A13299 Phthaldialdehyde, 98%

Product Number: A13299
CAS number: 643-79-8
MDL number: MFCD00003335
Molecular formula: C₈H₆O₂
Molecular weight: 134.14



Product Specification

Appearance (Color): Pale yellow to yellow
Form: Crystals, crystalline powder or lumps
Assay (GC): ≥97.5
Free acid (titration): ≤1.0% (as Phthalic acid)
Identification (FTIR): Conforms
Water Content (Karl Fischer Titration): ≤0.75%
Melting Point (clear melt): 51-59°C

Date of Print: November 21, 2022
Version: 6

Product Specifications are subject to amendment and may change over time.