



Product Specification Sheet

VWR MOP HEADS FOR PERFEX MOP SYSTEM

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No-US	Cat. No-EU	Description	Pk
82024-316	129-0047	PU foam	50
82024-318	129-0048	PU foam, tear-proof, extra absorbent	50

Other Features:

- Proprietary polyurethane foam compatible with disinfectants used for maintaining sterile environments.
- Laminated versions allow tear-free cleaning and disinfecting of areas that have rough or sharp surfaces.
- Micro-fiber versions are designed to remove and absorb stubborn contamination.
- Individually packaged to maintain cleanliness.
- Laminated versions are designed to allow good liquid transfer.
- Available with Lot Certification of Conformance for ease of audit and validation procedures.

Autoclave Settings:

- Exposure Temperature: 121°C
- Exposure Time: 30 Minutes
- Cycle Type: Gravity
- Cycle Time: 20 Minutes
- Pressure: 17.5 psig

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0



Product Specification Sheet

VWR MOP HEADS FOR PERFEX MOP SYSTEM

Physical properties:

- Density: 2.0 lbs./ft³ (Average)
- Pore Size: 90 PPI (Average)
- Tear Strength: 3.5 lbs./in (Average)
- Tensile Strength: 25 PSI (Average)

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use- If there is any doubt ask us.

Date of issue: **November 2012**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0