


Annexin V (APC) (A285884)

Specifications:

Name:	Annexin V (APC)
Applications:	Flow Cytometry
Dilutions:	Ready-to-use solution.
Expression System:	Escherichia coli
Nature:	Recombinant
Conjugate:	APC
Purity:	> 99% (by SDS-PAGE).
Purification:	Bacterially expressed Annexin V (APC) is purified to homogeneity and conjugated under optimum conditions. The conjugate is purified by size-exclusion chromatography.
Product Form:	Liquid
Formulation:	Supplied in Phosphate Buffered Saline, pH 7.4, with 15 mM Sodium Azide.
Storage:	Store in the dark at 2-8°C. Avoid prolonged exposure to light. Do not freeze!
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:



Separation of Annexin stained apoptotic Jurkat cells (red-filled) from unstained cells (black-dashed) in flow cytometry analysis (surface staining) of camptothecin induced apoptotic Jurkat cells stained using Annexin V (APC) (5 µl reagent per million cells in 100 µl of cell suspension).