



ZPACE8000477 PLASTIC DRUM LIQ 25 L DIN 61 UN



SPECIFICATIONS

UN-certification :	3H1/Y1.4/200
Material :	HDPE Black
Neck Type :	DIN 61
Height :	371,00 mm
Length :	384,00 mm
Wide :	289,00 mm



Pump It



Stainless steel pump

AVAILABLE ACCESSORIES

Key for opening :	29548.053
Tap :	29548.043
Stainless steel Pump :	BURK5601-0400
Stainless steel Pump:	BURK5601-0401
AccuOne Pump:	BURK5900-1050
PTFE Pump:	BURK5606-0600
Pump it	BURK5202-1000



PTFE pump



AccuOne pump

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343