

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Sample Dippers
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2024-09-30

Sample dipper, PTFE, 600 mm

Catalogue Number	Diameter Ø [mm]	Capacity [mL]	Pk
331-0006	54	100	1
331-0007	66	250	1
331-0008	80	500	1
331-0009	100	1000	1

Replacement shaft

Catalogue Number	Description	Material	Length [mm]	Capacity [mL]	Pk
331-0005	Replacement/extension shaft	PTFE	600	-	1

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Virgin Polytetrafluoroethylene (PTFE). Raw Material complies with Directive FDA 21 CFR 177.1550 and complies with USP Class VI - USP 27, NF 22,2004 for plastic at 70 °C.
Temperature Resistance	–200 °C to +280 °C.
Sterilization	Non-Sterile.
ATP Assay	Not tested.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not tested.
Non-Pyrogenic Statement	Not tested.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	DEHP are not used in the manufacture of the raw material and are not expected to be present.
RoHS	Raw materials does not contain any of the restricted substance above the maximum concentration values as listed in Annex II of RoHS 2 directive 2011/65/EU as amended by EU 2015/863.
REACH Statement	To the best of our knowledge material used for the production does not contain more than 0.1 % (w/w) of any substance of very high concern.

Storage Conditions Ambient Temperature.

Shelf Life Indefinite when stored under ideal conditions.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	2024-09