



Retort Stand Base, stainless steel

thickness 8 mm

ground surface

M 10 threaded

adjustable foot

rubber feet

Article no.	length mm	width mm	
2012	210	130	1700 g
2013	250	160	2500 g
2014	300	150	2800 g
2015	315	200	4000 g



© 2015

Carl Friedrich Usbeck KG
Industriestr. 12-14
42477 Radevormwald
Germany

Phone: +49 (02195) 91 18-0
Telefax: +49 (02195) 91 18-40

<https://www.usbeck.eu>