

Material	25654.361
Material description	Methyl ethyl ketone
Grade	TECHNICAL
CAS Number	78-93-3
Molecular formula	C ₂ H ₅ COCH ₃
Molecular mass	72.11

Characteristics	Specifications
Assay	≥ 99 %
Boiling point	79 - 80 °C
Density (20/4)	0.804 - 0.807

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.