

VWR® TRACEABLE® FLASHING LED BIG-DIGIT DUAL CHANNEL TIMER

Compact, multi-tasker
with 6-digit, dual channel
display

Individually calibrated

ISO 17025 certified

Individually certified—
traceable to NIST



Distinct alarms for both activities

- Features: two simultaneous countdown activities, flashing LED timing indicator/alarm, adjustable volume
- Unique visual indication for each channel; green light flashes while counting down and red light flashes once zero has been reached
- Adjustable volume: high decibel, low decibel, or off for silent flashing LED operation
- Large 0.625-inch high LCD digits may be read from 20 feet

SPECIFICATIONS

- Flashing LED timing indicator/alarm
- Adjustable volume
- Size: 3.5 x 3 x 0.625"; weight: 3 oz.
- Supplied: Traceable Certificate, magnetic back, flip-open stand/clip, opening for lanyard, and batteries

Range	Resolution	Accuracy	Channels	Cat. No.
99 hours, 59 minutes, 59 seconds	1 second	0.01%	2	89140-196

Additional VWR® Traceable® Timing Products

VWR® Traceable Three-Channel Alarm Timer

Innovative timer has programmable alarm with adjustable volume control anywhere

- Features: count up/down, memory recall, alarm, clock, adjustable alarm duration/volume
- Triple display shows three different times at once; tests can be started individually or simultaneously
- Adjustable volume alarm can be set to sound for 60 seconds or to sound continuously until cleared
- Supplied: battery, Traceable Certificate

Timing Capacity	Resolution	Accuracy	Channels	Cat. No.
99 hours, 59 minutes, 59 seconds	1 second	0.001%	3	62344-600



VWR® Traceable Vibrating Timer

Innovative countdown timer vibrates when time is up; quietly nudges you during a meeting or perfect for loud environments when hearing a buzzer is impossible

- Features: countdown/up, memory, stopwatch mode, water-Resistant, steamproof
- Individually numbered keys make setting time fast and easy
- Select either vibrating or an audible alarm with the flip of a switch
- Supplied: magnetic back, stand/spring-clip, Traceable Certificate, battery, and opening for lanyard

Timing Capacity	Resolution	Accuracy	Channels	Cat. No.
9 hours, 99 minutes, 99 seconds	1 second	0.01%	1	10154-394



VWR® Traceable Flashing LED Alert Big-Digit Alarm Timer

Big-Digit Timer has it all—adjustable alarm volume, flashing led, superb accuracy

- Features: count up/down, flashing LED, adjustable volume
- Green light flashes every second after timing begins; red LED indicator flashes when zero is reached and shows time elapsed since alarm
- Adjustable volume: high decibel, low decibel, or off for silent visual flashing LED operation
- Supplied: batteries, magnetic back, flip-open stand/ clip, opening for lanyard, Traceable Certificate

Timing Capacity	Resolution	Accuracy	Channels	Cat. No.
99 minutes, 59 seconds	1 second	0.01%	1	89140-194



Traceable to NIST— Multi-point calibration on an individually-numbered Traceable Certificate which assures accuracy from an ISO/IEC 17025:2005 (1750.01) calibration laboratory accredited by A2LA. It indicates traceability of measurements to the SI units through NIST or other recognized national measurement institutes (NMI) that are signatories to the CIPM Mutual Recognition Agreement.



1.800.932.5000 | vwr.com

Prices and product details are current when published; subject to change without notice. Certain products may be limited by federal, state, provincial, or local regulations. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.