



Product Specification Sheet

VWR GAS TIGHT MICROSYPINGES

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No: **549-0536**

Volume : **1 ml**

Measurement

Item	Standard
length from head to flange	98±2
barrel length	81.5±2
barrel O.D.	φ10±0.4
length from tip end to flange top	97±2
length of plunger ASSY	106.5±2
direction of flange printing	printed at the opposite of milled flange

Part	Material
PTFE Lure tip	PTFE
Lure head	BS+Cr coating
Barrel	hard glass
Flange	SUS303
Plunger	SUS304
Knob	SUS303
plunger tip	PTFE



Product Specification Sheet

VWR GAS TIGHT MICROSYRINGES

Packing method : **packed in paperboard box, urethane(inside)**

Storage condition :

Temperature : **from -10°C to 5**

To prevent change in resin part shape, deterioration of bond and crack at adhesive part by heat expansion.

Humidity : **from 20% to 90%**

Dust: **Keep away from the sticking of dust.**

Ultraviolet : **Keep away from UV.**

Warranty expiration : **36 months under the above mentioned condition**

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt, ask us.