

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Septa ND11
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2024-01-03

Catalogue Number	Septum	Hardness [shore A]	Thickness [mm]	Pk
548-0395A	Cream silicon/red PTFE	55°	1.5	1
548-0843A	Grey PTFE/cream butyl/grey PTFE	55°	1.3	1
548-1554A	White silicone/silver aluminium	50°	1.3	1

Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Visual Inspection: No deformation, contamination, enclosure.

Product Specifications

Material	Refer table.
Sterilization	Non-Sterile.
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.

BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	See RoHS.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Dry, protected from light and standard warehouse temperature.
Shelf Life	2 years, conditions: store clean, dry, and protected from light at normal storage temperature.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	10-2024