

Certificate of Quality

VWR NA Catalogue Number	VWRI210-000002
VWR EU Catalogue Number	VWRI210-000002
Description	VWR®, Cryovials, External threaded vials, Natural, 1.00/1.20 mL
Lot Number	VFAS13C25S
Date of Manufacture (yyyy-mm)	2025-03
Date of Sterilization (yyyy-mm)	2025-03
Sterilization Expiry Date (yyyy-mm)	2030-02
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2025-08-26

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 13485:2016**

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-03	General Insp. Level-1, AQL -1.5	Leak Proof
Volume Accuracy	QWI-QC-13	General Insp. Level-1, AQL -1.5	Within tolerance
Radiation Dose	QWI-QC-02	General Insp. Level-1, AQL -1.5	PASS (Minimum Dose- 15 K Gray)
Dimensional	QWI-QC-23	General Insp. Level-1, AQL -1.5	Within tolerance

Product Specifications

Products are manufactured in cleanroom environment (Class 10,000).

Material	100 % Virgin Tube Material: Polypropylene Cap material: High Density Polypropylene
Temperature Resistance	-196 °C.
Sterilization	Products are Gamma irradiated as per ISO 11137 and meets SAL 10 ⁻⁶ .
ATP Assay	Not Applicable.
DNase & RNase Free	The product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Products are free of cytotoxicity effects.
Non-Pyrogenic Statement	Products are free from pyrogen.
Latex Statement	The products are not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.

Storage Conditions	Storage and Transport Conditions-Store and transport at ambient Temperature. Do not expose to extreme temperature and direct Sunlight.
Shelf Life	5 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.